

Individual 800x macro stocks	Grams/L
NH ₄ NO ₃	91.4
NaH ₂ PO ₄ · 2H ₂ O	40.3
K ₂ SO ₄	71.4
CaCl ₂	88.6
MgSO ₄ · 7H ₂ O	324

Complete 800x micro stock	Grams/L
MnCl ₂ · 4H ₂ O	1.5
(NH ₄) ₆ MO ₇ O ₂₄ · 4H ₂ O	0.074
H ₃ BO ₃	0.934
ZnSO ₄ · 7H ₂ O	0.035
CuSO ₄ · 5H ₂ O	0.031
FeCl ₃ · 6H ₂ O	7.7
Citric acid (monohydrate)	11.9

micro stock instructions

dissolve separately:

- i) 1.5g MnCl₂·4H₂O in 40 mL H₂O bring to 50 mL
- ii) 0.074g (NH₄)₆MO₇O₂₄·4H₂O in 40 mL dH₂O bring to 50 mL
- iii) 0.934g H₃BO₃ in 40 mL bring to 50 mL
- iv) 7.7g FeCl₃·6H₂O in 40 mL H₂O bring to 50mL (do this in the hood)
- v) 11.9g Citric acid monohydrate in 40mL H₂O bring to 50 mL
- vi) prepare [50mM] stock of ZnSO₄ · 7H₂ O --> 14.37 g/ Liter
- vii) prepare [50mM] stock of CuSO₄ · 5H₂ O --> 12.48 g / Liter
- viii) add 50mL conc. H₂SO₄ to 550mL dH₂O in glass beaker
- ix) slowly add solutions to diluted H₂SO₄:
 - 1.) 50 mL of MnCl₂·4H₂O
 - 2.) 50 mL of (NH₄)₆MO₇O₂₄·4H₂O
 - 3.) 50 mL of H₃BO₃
 - 4.) 2.435 mL of [50mM] ZnSO₄·7H₂O
 - 5.) 2.483 mL of [50mM] CuSO₄·5H₂O
 - 6.) 50 mL of FeCl₃·6H₂O (do this in the hood)
 - 7.) 50 mL of Citric acid monohydrate
 - 8.) Bring up to 1L

Yoshida's gel for rice

1. add 1.25 mL each macro stock solution and the micro stock to 1 Liter final volume
2. add 0.546 g MES hydrate per Liter
3. pH to 5.8 with NaOH
4. Add Gelzan (Caisson) 2.5g/ Liter for a final 0.25% gel
5. Autoclave with stir bar and mix prior to gel hardening